


Innovative Reliable Quality
Manufacturer & Exporter of Surgical Instruments

- Stainless Steel
- Titanium



NAZIR & HABIB
Corporation

OPHTHALMOLOGY

Edition 2023-2024

Correspondence:

1285 Nakhawal, Kotli Loharan West Sialkot (51310)

Mob/WhatsApp: +92 346 6633 121, +92 317 7755121

Email: nazirhabibcorporation@gmail.com, sales@nazirhabibcorp.com

Web: www.nazirhabibcorp.com

Bank Alfalah Limited, Sialkot Pakistan (NTN-4238139-8)

Pakistan (51220)





NAZIR & HABIB
Corporation

OPHTHALMIC INSTRUMENTS

Innovative Reliable Quality



COMPANY PROFILE:

NAZIR & HABIB CORPORATION is manufacture and specialized in **MICRO OPHTHALMIC INSTRUMENTS**.

Using state of the art Technology and working in close collaboration with surgeon, **NAZIR & HABIB CORPORATION**

has designed, produced and micro finished the real micro ophthalmic surgeon instruments

Regards,

A.MANNAN
(C.E.O)





GENERAL INFORMATION

Warranty/Guarantee:

All NAZIR & HABIB Corporation instruments are un conditionally guaranteed for life against manufacturing defects in Material and workmanship when used for their intended surgical purpose and cared for in accordance with recommended procedure.

In the event any NHC instruments should fail as a result of a manufacturing defect in material or workmanship, It will be replaced or repair at no charge.

TERMS OF PAYMENTS:

As per invoice payment will be 100% advance for newly client and also note that after the one year relationship and regular orders our policy will be changed 75% advance and 25% upon shipment.

SHIPPING:

Shipment of surgical instruments are made via DHL, Fedex and EMS. Also note that upon the customer request, with a valid shipping account number for requested shipping vendor. If customer use their shipping account number a handling fee will still apply to the order.

ORDERS:

Orders can be placed through our emails: nazirhabibcorporation@gmail.com, sales@nazirhabibcorp.com and for the quick response you can send inquiry our whatsapp no: 0092 346 6633121, 0092 317 7755121.

PRICES:

Prices are in US Dollars and are subject to change without notice.

NOTE:

The following items CANNOT be accepted for credit excepts those resulting from an error made by NHC.

- Custom made instruments
- Discontinued models
- Damaged, abused, or etched instruments
- Sterile and disposable products
- Cannulas
- Diamond Knives



TABLE OF CONTENTS

TITANIUM INSTRUMENTS SECTIONS


▲ TITANIUM FORCEPS	PAGE 1
▲ TITANIUM SCISSORS	PAGE 2
▲ TITANIUM NEEDLE HOLDERS	PAGE 3
▲ TITANIUM HOOKS	PAGE 4-5
▲ TITANIUM CLAMPS	PAGE 6
▲ TITANIUM CANNULAS	PAGE 7

STAINLESS STEEL INSTRUMENTS SECTIONS

▲ STAINLESS STEEL SCISSORS	PAGE 8
▲ STAINLESS STEEL FORCEPS	PAGE 9-10
▲ STAINLESS STEEL NEEDLE HOLDERS	PAGE 11
▲ STAINLESS STEEL HOOKS	PAGE 12-13-14
▲ STAINLESS STEEL CLAMPS	PAGE 15
▲ STAINLESS STEEL SPECULUMS	PAGE 16-17
▲ STAINLESSS STEEL DCR	PAGE 18-19
▲ MISCELLANEOUS	PAGE 20
▲ IRRIGATION/ASPIRATION	PAGE 21
▲ CANNULAS & PROBES	PAGE 22 -23
▲ CORNEAL MARKERS	PAGE 24
▲ VR INSTRUMENTS RANGE	PAGE 25-26-27-28-29-30-31
▲ HOLLOWWARE	PAGE 32




FORCEPS TITANIUM



NHC: 786-0201



Corneal Forceps Straight, Titanium
1x2 teeth, 0.1 mm height, Tips 0.3 mm wide



NHC: 786-0202



Castroviejo Tying Forceps Straight, Titanium
1x2 teeth 0.1 mm height, 6.0 mm tying platforms



NHC: 786-0203




Utrata Capsulorhexis Forceps, Angled, Titanium
Overall length 102 mm, 4 inches, Delicate shafts 10 mm long



NHC: 786-0204



Masket Capsulorhexis Forceps very delicate
triangular grasping tips extremely thin
11mm long vaulted shanks



NHC: 786-0205



McPherson Forceps Angled 7.5mm long jaws



NHC: 786-0206



McPherson Tissue Forceps 1x2 teeth, 0.4mm
with tying platform



NHC: 786-0306




Colibri Forceps standard 1x2 teeth 0.4mm




NHC: 786-0207




KLW
Bonn Forceps very delicate 1x2 teeth 0.12mm




NHC: 786-0208



Dressing Forceps, Serrated delicate straight




NHC: 786-0209



Tissue Forceps 1x2 teeth delicate curved



NHC: 786-0210
Lieberman Eye Speculum



NHC: 786-0211
Lieberman Tanner Eye Speculum

Jewelers #5 Forceps, 11cm, Straight, Titanium
NHC: 786-0212



NHC: 786-0213
Jewelers #7 Forceps, 11cm, 0.07x0.01mm Tips, Titanium



SCISSORS TITANIUM



NHC: 786-0214
Corneal Scissor Curved, Blunt tips



NHC: 786-0215
Westcott Tenotomy, Scissors Curved blunt tips



NHC: 786-0216
Stitch Scissor, Sharp blade, Straight




NHC: 786-0217
Stitch Scissor, Sharp blade, Curved



NHC: 786-0218
Corneoscleral Scissors, Medium blades Right



NHC: 786-0219
Corneoscleral Scissors, Medium blades Left



NHC: 786-0220
Vannas Scissors, Sharp tips Straight



NHC: 786-0221
Vannas Scissors, Sharp tips Curved




NHC: 786-0222
Vannas sharp tips, extra thin angled forward




NHC: 786-0223
Vannas Sharp tips, angled to side




NHC: 786-0224
Vannas Capsulotomy Scissors, 9.5cm, 10.5mm
Sharp Angled Tips



NHC: 786-0225
Castroviejo Universal Corneal Scissors
curved, blunt tips, small blades




NHC: 786-0226
Stevens Tenotomy Scissors
blunt tips, standard Straight




NHC: 786-0227
Eye Scissor Standard, Pointed tip
Curved and Straight available




NEEDLE HOLDER TITANIUM




NHC: 786-0228
Needle holder delicate jaws curved without lock



NHC: 786-0229
Needle holder micro curved without lock



NHC: 786-0230
Needle holder micro extra delicate jaws straight



NHC: 786-0231
Needle holder micro curved with lock




NHC: 786-0232
Castroviejo needle holder curved without lock




NHC: 786-0233
Castroviejo needle holder straight without lock




NHC: 786-0234
Castroviejo needle holder curved with lock



NHC: 786-0235
Castroviejo needle holder straight with lock




NHC: 786-0236
Barraquer Micro needle holder short model




NHC: 786-0237
Troutman Needle Holder
delicate jaws, strongly curved without lock




NHC: 786-0238
Micro Castroviejo Needle Holder, 9cm, Straight



NHC: 786-0239
Micro Needle Holder Flat Handle



NHC: 786-0240
Hartmann mosquito forceps extra
delicate Curved



NHC: 786-0241
Halsted mosquito forceps delicate Straight



HOOKS TITANIUM



NHC: 786-0242

Sinskey lens manipulating hook with round handle 0.2mm dia. Blunt tip angled



NHC: 786-0243

Lewicky lens manipulating hook with round handle 0.15mm dia. Blunt tip angled



NHC: 786-0244

Bechert Lens pusher vertical Y shaped tip angled



NHC: 786-0245

Alio-Kelman
IOL Hook 0.2mm diameter tip



NHC: 786-0246

Iris Spatula
angled 2mm wide x 8mm long 0.2mm thick



NHC: 786-0247

Osher Y-Hook lens Manipulator straight with guard



NHC: 786-0248

Milder Phaco Spatula
60° offset vertical and horizontal chopper



HOOKS TITANIUM



NHC: 786-0249
Nagahara Phaco Chopper rounded blunt tip 0.9mm inferior wedge



NHC: 786-0250
Phaco Chopper Ball tips 0.7mm inferior wedge



NHC: 786-0251
Minardi Phaco Chopper 0.75mm pointed "quick chop" tips



NHC: 786-0252
Hunkeler IOL Manipulator with 0.3mm diameter hemi ball tip



NHC: 786-0253
Fenzl Lens Manipulating hook 0.15mm diameter tips




NHC: 786-0254
Kuglen Iris hook and manipulator push pull model round handle




NHC: 786-0255
Lens Loop




CLAMPS TITANIUM



NHC: 786-0256
Lambert Chalazion Forceps Round 8mm



NHC: 786-0257
Baird Chalazion Forceps oval 11 x 8mm I.D small



NHC: 786-0258
Francis Chalazion Clamps oval, Medium



NHC: 786-0259
Akahoshi combo prechopper for cracking
dividing and rotating



NHC: 786-0260
Titanium Vitreoretinal Forceps
23ga Straight Blunt Tips



NHC: 786-0261
Capsulorrhexis Forceps
delicate tapered shanks



NHC: 786-0262
Inamura
Capsulorrhexis Forceps
Direct Action style 3mm
incision size






CANNULAS TITANIUM


NHC: 786-0263

Simcoe Cannula 1 Irrigating Aspirating Cannula
15mm long, 0.3mm diameter, 23 gauge thin wall




NHC: 786-0264

Reversed Simcoe Cannula 1 Irrigating Aspirating Cannula
15mm long, 0.3mm diameter, 23 gauge thin wall



NHC: 786-0265

Gess Irrigating Aspirating Cannula U-shaped aspiration tip, end opening port 23 gauge




NHC: 786-0266

Air Injection Cannula angled at 5mm



NHC: 786-0267

Bowman Lacrimal Probe Malleable Stainless Steel




NHC: 786-0268

Wilder Lacrimal Dilator





SCISSORS




NHC-10100
Corneal Scissor Curved, Blunt tips



NHC-10101
Westcott Tenotomy, Scissors Curved blunt tips



NHC-10102
Stitch Scissor, Sharp blade, Straight



NHC-10103
Stitch Scissor, Sharp blade, Curved



NHC-10104
Corneoscleral Scissors, Medium blades Right



NHC-10105
Corneoscleral Scissors, Medium blades Left



NHC-10106
Vannas Scissors, Sharp tips Straight



NHC-10107
Vannas Scissors, Sharp tips Curved




NHC-10108
Vannas sharp tips, extra thin angled forward



NHC-10109
Vannas Sharp tips, angled to side




NHC-10110
Keratoplasty Scissors angled to side, blunt tips, medium blades



NHC-10111
Eye Scissors Standard, Pointed tip



NHC-10112
Enucleation Scissors medium and strong curve blunt tips



STAINLESS STEEL



FORCEPS




NHC-10113
Tying Forceps, 5mm long smooth jaws straight



NHC-10115
Corneal Forceps 1x2 teeth with tying platform



NHC-10117
McPherson Forceps Angled 7.5mm long jaws



NHC-10119
Utrata Sapsulorhexis Forceps 11mm long shanks



NHC-10121
Tissue Forceps 1x2 teeth delicate straight



NHC-10123
Dressing Forceps, Serrated delicate straight



NHC-10125
Arruga Capsule Forceps



NHC-10114
Tying Forceps, Standard Smooth jaws curved



NHC-10116
Bonn Forceps very delicate 1x2 teeth 0.12mm



NHC-10118
Colibri Forceps standard 1x2 teeth 0.4mm



NHC-10120
Capsulorhexis Forceps 11mm Vaulted shanks



NHC-10122
Tissue Forceps 1x2 teeth delicate curved



NHC-10124
Dressing Forceps, Serrated delicate curved



NHC-10126
Utility Forceps very delicate smooth jaws curved

STAINLESS STEEL



FORCEPS



NHC-10127
Akahoshi combo prechopper for cracking
dividing and rotating



NHC-10128
Nucleus cracker cross action



NHC-10129
Clayman lens holding forceps



NHC-10130
Clayman lens holding forceps



NHC-10131
Lens Holding Forceps 4.5" angled round smooth
polished jaws.



NHC-10132
Lens Folding Forceps cross action



NHC-10134
Humer Lens Folding Forceps



NHC-10133
Lens Folding Forceps with mid jaw indicators



NHC-10135
Sleeve Spreading Forceps



NHC-10136
Kelly Descemets Membrane Punch



NHC-10137
Trabeculectomy punch angled cutting blade



NHC-10138
Micro Ring Lens Forceps 1.2mm ring



NHC-10139
Bipolar Forcep Straight




NHC-10140
Bipolar Forcep Curved

STAINLESS STEEL



NEEDLE HOLDERS



NHC-10141
Needle holder delicate jaws curved without lock



NHC-10143
Needle holder micro extra delicate jaws straight



NHC-10145
Castroviejo needle holder curved without lock




NHC-10147
Castroviejo needle holder curved with lock




NHC-10149
Barraquer Micro needle holder short model



NHC-10151
Swiss model blade holder pencil type




NHC-10153
Hartmann mosquito forceps extra delicate



NHC-10142
Needle holder micro curved without lock



NHC-10144
Needle holder micro curved with lock



NHC-10146
Castroviejo needle holder straight without lock




NHC-10148
Castroviejo needle holder straight with lock



NHC-10150
Serrefine Serrated Jaws small straight



NHC-10152
Silcock needle holder



NHC-10154
Halsted mosquito forceps delicate

STAINLESS STEEL



HOOKS



NHC-10155

Sinskey lens manipulating hook with round handle 0.2mm dia. Blunt tip angled



NHC-10156

Lewicky lens manipulating hook with round handle 0.15mm dia. Blunt tip angled



NHC-10157

Bechert Lens pusher vertical Y shaped tip angled



NHC-10158

Crewford Retriving Hook



NHC-10159

Helveston Muscle Hook 8mm, 10mm & 12mm
extra delicate, with blunt spatulated tip



NHC-10160

Osher Y-Hook lens Manipulator straight with guard



NHC-10161

Milder Phaco Spatula
60° offset vertical and horizontal chopper

STAINLESS STEEL



HOOKS



NHC-10162
Nagahara Phaco Chopper rounded blunt tip 0.9mm inferior wedge



NHC-10163
Phaco Chopper Ball tips 0.7mm inferior wedge



NHC-10164
Minardi Phaco Chopper 0.75mm pointed "quick chop" tips



NHC-10165
Hunkeler IOL Manipulator with 0.3mm diameter hemi ball tip



NHC-10166
Fenzl Lens Manipulating hook 0.15mm diameter tips



NHC-10167
Kuglen Iris hook and manipulator push pull model round handle




NHC-10168
Drysdale Nucleus Manipulator
paddle shaped tip

STAINLESS STEEL



HOOKS

NHC-10169
Bunge Evisceration Spoon




NHC-10170
Wells Enucleation Spoon



NHC-10171
Fixation Hook Double, Sharp



NHC-10172
Graefe Muscle Hook size 1-2 & 3



NHC-10173
Gass Retinal Detachment Hook with oval hole



NHC-10174
Jameson Muscle Hook Small & Large



NHC-10175
Braunstein Fixed Caliper
3.5mm - 4.0mm



NHC-10176
Caliper

STAINLESS STEEL



CLAMPS



NHC-10177
Hunt Chalazion Forceps round, 12mm I.D.




NHC-10178
Lambert Chalazion Forceps Round 8mm



NHC-10179
Wies Chalazion Forceps grooved upper jaw



NHC-10180
Baird Chalazion Forceps oval 11 x 8mm I.D small



NHC-10181
Desmarres Chalazion Forceps
oval, 20 x 13mm I.D. large



NHC-10182
Francis Chalazion Clamps oval, Medium



NHC-10183
Snellen Entropium Clamp left



NHC-10184
Snellen Entropium Clamp right

NHC-10185
Berke Ptosis Clamps with side lock longitudinally grooved jaws



NHC-10186
Backhaus Towel Clamp




NHC-10187
Bobby Jones Towel Clamp cross-action




STAINLESS STEEL



SPECULUMS




NHC-10188
Kratz Barraquer Wire Speculum with open blades




NHC-10190
Barraquer Wire Speculum with solid blades




NHC-10192
Infant Wire Speculum




NHC-10194
Infant Wire Speculum



NHC-10189
Barraquer Wire Speculum



NHC-10191
Barraquer Wire Speculum



NHC-10193
McIntyre Wire Speculum
gently curved springs




NHC-10195
Alfonso Infant Wire Speculum

STAINLESS STEEL




SPECULUMS




NHC-10196
Lieberman Eye Speculum



NHC-10197
Lieberman Tanner Eye Speculum




NHC-10198
Murdoch eye Speculum fenestrated 9mm blades.



NHC-10199
Fishkind Eye Speculum
open blades adjustable mechanism



NHC-10200
Lieberman eye speculum



NHC-10201
Lancaster eye speculum solid blades



NHC-10202
Castroviejo Speculum Open Blades



NHC-10203
Sauer Lid Speculum for infant 11mm blades

STAINLESS STEEL



DCR

NHC-10204
Kerrison Bone Nibbling Rongeur, 9mm Opening



NHC-10205
Beyer Rongeur, single action



NHC-10206
Nasal Speculum



NHC-10207
Bone Rongeur, Double action



STAINLESS STEEL



DCR



NHC-10208

Freer Periosteal Elevator one end blunt other end semi-sharp 4.5 mm wide blades



NHC-10209

West Lacrimal Chisel Straight



NHC-10210

Pigtail Probe with suture holes



NHC-10211

Mallet filled stainless steel head



NHC-10212

Knapp Lacrimal Retractor 4 prongs



NHC-10213

Stevenson Lacrimal Sac Retractor 3 x 3 blunt prongs, adjustable



NHC-10214

Bard Parker Handle
#3, for disposable slot blades



STAINLESS STEEL



MISCELLANEOUS



NHC-10215
Schocket Scleral Depressor double-end with pocket clip




NHC-10216
Desmarres Lid Retractor thin solid blades



NHC-10217
Jaegar Lid Plate stainless steel



NHC-10218
Wright Fascia Needle 1mm x 6mm oval hole 140mm long



NHC-10219
Meyerhoefer Chalazion Curette



NHC-10220
Wheeler Double-End Spatula lightly curved



NHC-10221
Lens Loop

STAINLESS STEEL



IRRIGATION-ASPIRATION HAND PIECES



NHC-10222

Irrigation/Aspiration 45 Degree Handpiece

1-Piece Design With An Overall Length Of 5 1/2" (14cm). Overall Weight of 1.8 oz.



NHC-10223

Irrigation Aspiration system titanium 20g, 21g, 23g

Lightweight, color-coded irrigation and aspiration handpieces

Compatible with your existing Phaco or I/A System



NHC-10224

Irrigation Aspiration system stainless steel 20g, 21g, 23g

Lightweight, color-coded irrigation and aspiration handpieces


Compatible with your existing Phaco or I/A System



NHC-10225

Alio Irrigating Stinger

oval, 20 gauge



NHC-10226

Strauss-Simcoe

**Irrigation/Aspiration Cannula 23 gauge, 0.3mm port
elongated and gently curved with 15° angled tip**


STAINLESS STEEL



CANNULAS & PROBES


NHC-10227

Simcoe Cannula 1 Irrigating Aspirating Cannula
15mm long, 0.3mm diameter, 23 gauge thin wall




NHC-10228

Reversed Simcoe Cannula 1 Irrigating Aspirating Cannula
15mm long, 0.3mm diameter, 23 gauge thin wall




NHC-10229

Gess Irrigating Aspirating Cannula U-shaped
aspiration tip, end opening port 23 gauge




NHC-10230

Air Injection Cannula angled at 5mm




NHC-10231

Bowman Lacrimal Probe Malleable Stainless Steel



NHC-10232

Wilder Lacrimal Dilator





STAINLESS STEEL



CANNULAS & PROBES

NHC-10233
Lacrimal Cannula
malleable tip 23 gauge



NHC-10234
Troutman Olive Tip Cannula
delicate 25 gauge



NHC-10235
Vidaurri LASIK Cannula double-armed
8 irrigating ports 25 gauge



NHC-10236
Randolph Cyclodialysis Cannula
formed flat tip angled at 13mm



NHC-10237
Lacrimal Intubation Set
with retrieving hook
sterile, disposable



NHC-10238
Bodian
Lacrimal Probe Set
includes 2 curved blunt 1.0mm dia. probes
attached to 30.5cm long silicone tubing (0.012" I.D., 0.025" O.D.)



STAINLESS STEEL



CORNEAL MARKERS



NHC-10239
Mendez Degree Gauge
12mm I.D.



NHC-10240
Bores Axis Marker used with Mendez



NHC-10241
Toric IOL Marker
two 2.5mm radial marks extending to 11mm



NHC-10242
Greene Visual Axis Marker
0.5mm diameter concave tip



NHC-10243
Bores Optical Zone Marker
3.0mm, 3.25mm, 3.5mm, 3.75mm



NHC-10244
Anis Radial Marker
4.0mm center opening
12.5mm O.D. 8 blades, 12 blades



NHC-10245
Capsulorrhexis Marker 5.5mm diameter
marking ring with four
non-marking alignment peg



STAINLESS STEEL



VITREORETINAL RANGE

- Lightweight titanium handle reduces hand fatigue while providing lasting durability.
- Squeeze style handle provides extremely smooth actuation for high precision surgical manoeuvres.
- Rounded grip area enables easy 360° rotation during use for optimal positioning of tip.
- Forceps and scissors have stationary distal tips to eliminate any unintended movement of the jaws during actuation.
- Various sizes available including: 0.90mm (20 Gauge) 0.65mm (23 Gauge) 0.50mm (25 Gauge)
- A wide variety of instrument tips are available to meet every surgical need. Tips are made of unique Stainless Steel which is highly corrosion resistant, resulting in smoother operation and fewer repairs. Some tips are also available in Titanium for reduced glare and added




- Stainless Steel
- Titanium

20 GAUGE MICRO SCISSORS

VR-1010

90° Angled Scissors (Vertical)


Sharp tipped vertical style scissors with 1.9mm long blades.



VR-1011

Straight Scissors

Sharp tipped straight scissors with 3mm long blades.



VR-1012

Curved Scissors


Sharp tipped 3mm long blades curved to a 4mm radius.



VR-1013

135° Angled Scissors


Sharp tipped horizontal style scissors angled 135°, 2mm from the tip.



VR-1014

45° Angled Scissors (Vertical)


Sharp tipped vertical style scissors with 1.75mm long blades.



VR-1024

125° Angled Scissors


Sharp tipped horizontal style scissors angled at 125°, 2mm from the tip.



VR-1025

60° Angled Scissors (Vertical)


Sharp tipped fine vertical style scissors with 1.9mm long blades.



VR-1026

45° Scissors with Membrane Pick (Vertical)

Sharp tipped vertical style scissors with 2mm distal blade and 1.75mm proximal blade.



VR-1034

Vertical Curved Scissor

3mm long blades.





VITREORETINAL RANGE

SUB RETINAL AND ILM INSTRUMENTS

The treatment of Macular diseases can be complex and with new surgical approaches to macular conditions, this has meant that there is a need for a range of instruments which are more than capable of addressing these issues. Our ability to listen to surgeons, take ideas and then design an instrument means that we have a range of ILM and Subretinal instruments which are amongst the most advanced in the market.


VR-1035
Long Angled Scissors

Angled 135°, 3.5mm from the tip.




VR-1036
Angled Long Gripping Forceps

135°, 3mm from the tip.




VR-1037
Fine Gripping Forceps

3.5mm long jaws.



VR-1038
Ducournau Fine Gripping Forceps

Angled 120°, 3.5mm from the tip.
Developed in collaboration with Didier Ducournau,
M.D., Nantes, France.



VR-1050
Naito ILM Forceps

3.3mm long jaws with a 0.7mm pick –
Non Glare finish.



VR-1350
Naito ILM Forceps
(Titanium Jaws)

Developed in collaboration with
Dr. Naito, Assistant Professor,
Tokushima University
Hospital, Japan.




VR-1051
Naito Subretinal Forceps

Angled 135°, 4.5mm from the tip –
Non Glare finish




VR-1051s
Naito Subretinal Forceps

Angled 135°, 3mm from the tip – Non Glare finish.
Developed in collaboration with Dr. Naito, Assistant
Professor, Tokushima University Hospital, Japan.



VR-1053
Delicate Grasping Forceps

3.5 long, straight forceps.
Developed in collaboration with Professor
Jean-Paul Berrod, Centre Hospitalier
Universitaire, Nancy, France.





VITREORETINAL RANGE

SMALL INCISION TAPERED INSTRUMENTS VR-1700 SERIES

The VR-1700 series instruments are manufactured with a Tapered Shaft. These are designed with 20G proximal end with a distal end that tapers to 23G. The tips are some of the smallest instruments available for ophthalmic microsurgery.



VR-1703

Pointed End Gripping

Tapered Shaft (20G tapering to 23G).



VR-1712

Curved Scissors

Tapered Shaft (20G tapering to 23G)
2.25mm long blades curved to a
2.5mm radius.



VR-1711

Straight Scissors

Tapered Shaft - 2.25mm long blades
(20G tapering to 23G).



VR-1713

135° Angled Scissors


Tapered Shaft (20G tapering to 23G)
angled 135°, 1.5mm from the tip.






VITREORETINAL RANGE


VR-1530
23G Pick Forceps
3mm long jaws



VR-1539
23G Serrated ILM Forceps, 1.25mm
Long – with fine pick.
Designed for Fine Membrane and
manipulation of retinal folds.




VR-1531
23G Offset Gripping Forceps
0.9mm long jaws




VR-1550
23G ILM Forceps
Non Glare finish, 2.2mm long with shorter
shaft – 27.5mm long.



VR-1537
23G 3mm Long Fine Gripping
Forceps



VR-1553
23G Delicate Grasping Forceps
2.75mm long jaws






VITREORETINAL RANGE

ILLUMINATED MICRO INSTRUMENTS

Fibreoptic illuminated instruments give you the opportunity to improve visualisation by providing light at the tip. The special squeeze handle has the option of adding a 0.5mm flexible fibreoptic light pipe (available separately). The instruments utilize a 1.1mm (19 Gauge) shaft, which tapers to 0.65mm (23 Gauge) 3mm from the tip. This provides light where you need it most while freeing your other hand to hold a second instrument. The instrument tips are crafted in a unique, highly corrosion



VR-1112

Curved Scissors (Illuminated)

VR-1112

Curved Scissors (Illuminated)

2.25mm long blades curved to a 2.5mm radius.



VR-1118

Angled Forceps (Illuminated)

Angled 135°, 1.25mm from the tip with small platforms.



VR-1113

135° Angled Scissors (Illuminated)


Angled 135°, 1.5mm from the tip.



VR-1120

Serrated Forceps (Illuminated)

1.5mm long jaws.






VITREORETINAL RANGE

25G INSTRUMENTS


VR-1811
25G Straight Scissor
1.25mm long blade.




VR-1812
25G Curved Scissor
1.25mm long blade.




VR-1813
25G Long Curved Scissor
2.5mm long blade.




VR-1814
25G Vertical Scissor
0.5mm long blade.




VR-1816
25G Micro End Gripping
Forceps
1.25mm long jaws.




VR-1820
25G Serrated Forceps
1.25mm long jaws.




VR-1830
25G Pick Forceps
1.5mm long jaws.




VR-1831
25G Offset (Asymmetrical)
Gripping Forceps, 0.9mm
long jaws.



VR-1837
25G Fine Gripping Forceps
1.25mm long jaws.



VR-1850
25G ILM Forceps
Non Glare finish, with
shorter shaft –
27.5mm long.





VITREORETINAL RANGE

MICRO SCISSORS AND MICRO CAPSULORHEXIS FORCEPS

These micro instruments provide a new level of surgeon control. The Titanium squeeze handle provides delicate, tactile feedback and enables excellent tear control for Capsulorhexis. The small diameter shaft allows use through a side port or Phaco incision without stretching the incision or flattening the anterior chamber. The curved shaft allows easy access and comfortable hand position, even in deep set eyes.




MICRO CAPSULORHEXIS FORCEPS

FR-2262

Golenvaux Universal Capsulorhexis Forceps - 0.7mm long jaws with 90° angled grasping platform (23G)

Strongly curved shaft 18mm long.


(Developed in collaboration with Dr Golenvaux., Belgium)



FR-2266

Ikeda Capsulorhexis forceps 0.65mm (23G)

Delicate 1.2mm long jaws with angled, grasping platform. Curved shaft 24mm long (Developed in collaboration with Dr. Ikeda M.D., Japan).




FR-2265

Capsulorhexis forceps 0.80mm (21G)

4mm long jaws with angled, grasping platform.


Curved shaft 24mm long.



FR-2268

Ikeda angled jaws micro Capsulorhexis forceps 0.65mm (23G)

0.9mm long jaw with 90° angled, grasping platform. Strongly curved shaft 24mm long (Developed in collaboration with Dr. Ikeda M.D., Japan).




MICRO SCISSORS

FR-2275

2 mm long blades, 135 degree scissors 0.65mm (23G)

Small 0.65mm diameter, 32mm straight shaft allows access through a side port incision.



FR-2275C

2.5 mm long blades in curved shaft 0.65mm (23G)


Small 0.65 mm diameter, 24mm curved shaft allows access through a side port incision.






HOLLOWWARE

NHC-10246



Sterilization Box (Teflon Autoclavable 6 x 10 size)

NHC-10247



Sterilization Box (Teflon Autoclavable 4 x 7 size)

NHC-10248



Sterilization Box Stainless Steel

NHC-10249



Sterilization Cassettes

NHC-10250



Dressing Tray Stainless Steel

NHC-10251



Sterilization Dressing Drum

NHC-10252



Kidney Tray Stainless Steel

NHC-10253



Stainless Steel Bowls

Order Form (Please fill out completely)

Your Address

BILL TO

Attention: Accounts Payable

Phone ()

Order Number

SHIP TO

☐ Check if Same As Bill To Address

Signature

[illegible]

DHL 

FEDEX 

EMS ☐

407 Kotli Loharan West Sialkot (51310)

Mob: +92 346 6633121, +92 300 6133175

Email: salesklwinstruments@gmail.com, info@klwinstruments.com



NAZIR & HABIB Corporation

Dear Customer, Partner and our Friends,

NAZIR & HABIB Corporation is one of the leading manufacturers of high precision instruments and consumables of Ophthalmic surgery.

The company was established in 2012 we proud to be well known Surgical Instruments Company.

We are committed to provide surgeons with innovative and high quality instruments at reasonable prices. Our team is concerned about your patient's health that is why we are always follow the latest trends and innovative techniques in Surgeries.



For quick Inquiry:

Whatsapp no: +971 525 345 747, +92 332 62 90232

Email: nazirhabibcorporation@gmail.com